

Manual

Accompanying the GPHF-Minilab®

Supplement 2012

Volume II

THIN LAYER CHROMATOGRAPHIC TESTS



A charitable organisation
maintained exclusively by
Merck Darmstadt · Germany



USAID
FROM THE AMERICAN PEOPLE



PROMOTING THE QUALITY OF MEDICINES

A Concise Quality Control Guide on Essential Drugs and other Medicines

SUPPLEMENT 2012 TO VOLUME II ON THIN LAYER CHROMATOGRAPHIC TESTS

Written by

Richard W. O. Jähnke, Kornelia Dwornik, Kai Fischer and Souly Phanouvong

* * *

Reviewed by

Adrian Bajoras, Daniel Bempong, Sanford Bradby, Kennedy Chibwe, Yanga Dijiba, Latifa El Hadri, Mustapha Hajjou, Lukas Roth and Patrick Lukulay

* * *

Published by

the Global Pharma Health Fund (GPHF), a charity initiated and maintained by Merck Darmstadt · Germany, and the U.S. Pharmacopeia's Promotion the Quality of Medicines programme (USP PQM)

* * *

Copyright © by GPHF & USP PQM

* * *

Acknowledgements

This supplement was made possible by the generous financial support of the American people through the United States Agency for International Development (USAID). The contents are the responsibility of GPHF and USP PQM and do not necessarily reflect the views of USAID or the United States Government.

* * *

About the GPHF-Minilab® Project

Counterfeit medicines proliferation constitutes serious health hazards. The international police organisation Interpol estimates that a disturbing proportion of ten to thirty percent of all drugs offered in developing countries are either counterfeit or of deficient quality already. Fighting falsified medicines will ensure that decades of investments in healthcare are not undone through lack of vigilance.

To prevent counterfeit and extreme poor anti-infective medicines infiltrating drug supply organisations and priority disease programmes in malaria, TB and HIV/AIDS endemic countries, the Global Pharma Health Fund (GPHF) in Frankfurt, a charity maintained exclusively by Merck Darmstadt · Germany, set out to develop and supply at low cost the GPHF-Minilab®, a mini-laboratory for rapid drug quality verification and counterfeit medicines detection.

Since many years, GPHF-Minilabs are acting as a first-line defence against counterfeit and substandard quality medicines threatening the health of millions of people living in developing nations. Overall, more than 500 Minilabs have been supplied to over 80 countries across the African, Asian-Pacific and Latin American region already.

Main implementation partners are national health and medicines regulatory authorities together with the World Health Organization and the U.S. Pharmacopeia's Promoting the Quality of Medicines programme. Joint drug quality monitoring projects run by Interpol in South East Asia and East Africa triggered off the seizure of millions of counterfeit antimalarial pills without any active principles in the recent years.

The unchanged need for non-sophisticated and affordable drug quality monitoring in low-income countries forms the driving force behind the development of new GPHF-Minilab® test protocols today. The need for more testing emphasises the important collaboration with our US based implementing partners. For more patient safety and better health in developing countries, other parties are invited to join in.

* * *

Made by Grimm Graphic Design, Ochsenfurt, Germany

GPHF-Minilab® assembled and supplied by Technologie Transfer Marburg, Cölbe, Germany

Table of Contents

Chapter	Page
New TLC Test Protocols	4
<i>Supplement to Volume II, Chapter 6</i>	
<i>Artesunate-based fixed-dose combination antimalarial medicines:</i>	
6.58 Artesunate (incl. common co-formulations)	4
6.59 Amodiaquine (incl. artesunate co-formulations)	8
6.60 Mefloquine (incl. artesunate co-formulations)	12
6.61 Pyrimethamine (incl. SP and artesunate co-formulations)	16
6.62 Pyronaridine (as tetraphosphate incl. artesunate co-formulations)	20
6.63 Sulfadoxine (incl. pyrimethamine and artesunate co-formulations)	24
6.64 Sulfamethoxyprazine (incl. pyrimethamine and artesunate co-formulations)	28
Summary Table of Chromatographic Working Conditions	32
<i>Supplement to Volume II, Chapter 7</i>	
Updated List of GPHF-Minilab® Reference Standards	33
<i>Supplement to Volume II, Chapter 10</i>	
Health & Safety	35

6.58 Artesunate (incl. common co-formulations)

Primary Screening via Physical Inspection and Disintegration Test

I. PHYSICAL INSPECTION

Search for deficiencies on labelling, packaging and dosage forms as described in the opening chapters on general methods and operations of the main manual. Write down all product particulars using the reporting form as a guide. Having first been administered as a single drug, artesunate-based fixed-dose combination therapy is much more common now. Appropriate tablets or capsules are usually consisting of artesunate on one side and amodiaquine, mefloquine, pyronaridine, pyrimethamine, sulfadoxine or sulfamethoxyprazine on the other; the artesunate fraction usually being presented in dosage strengths between 25 and 200 mg.

II. DISINTEGRATION TEST

All quick release artesunate single oral dosage forms and appropriate fixed-dose combination tablets and capsules must pass the disintegration test as described in the opening chapters on general methods and operations of the main manual. They should disintegrate in water at 37 °C in less than 30 minutes. It's a major defect if a drug product doesn't pass this test.

III. RESULTS & ACTIONS TO BE TAKEN

Drug products from unusually cheap sources, drug products with missing or incorrect accompanying documents and drug products with defective dosage forms, packaging or with incomplete, damaged or missing labels or with labels written in a foreign language should be subjected to a thin layer chromatographic test.

Verification of Drug Identity and Content via Thin Layer Chromatography

I. PRINCIPLE

Whether or not combined with other drugs, artesunate is extracted from tablets and capsules with methanol and determined by TLC with reference to an authentic secondary standard. All procedures presented in this protocol are fit for the detection of artesunate in single oral dosage forms and appropriate fixed-dose combination products. For rapid drug quality verification of the amodiaquine, mefloquine, pyronaridine, pyrimethamine, sulfadoxine or sulfamethoxyprazine fraction consult the other protocols shown in this supplement.

II. EQUIPMENT AND REAGENTS

- 1) Pestle
- 2) Aluminium foil
- 3) Funnel
- 4) Label tape
- 5) Marker pen
- 6) Pencil and ruler
- 7) 10-ml vials
- 8) Set of straight pipettes (1 to 25 ml)
- 9) Set of laboratory glass bottles (25 to 100 ml)
- 10) Merck TLC aluminium plates precoated with silica gel 60 F₂₅₄, size 5x10 cm
- 11) Glass microcapillaries (2-µl filling capacity)
- 12) TLC developing chamber (500-ml jar)
- 13) Hot plate
- 14) Filter paper
- 15) Pair of scissors
- 16) Pair of tweezers
- 17) UV light of 254 nm
- 18) TLC dipping chamber (250-ml beaker)
- 19) Sulphuric acid solution 96%
- 20) Acetone
- 21) Ethyl acetate
- 22) Glacial acetic acid
- 23) Methanol
- 24) Secondary reference standard, for example artesunate 50 mg tablets

III. PREPARATION OF THE STOCK STANDARD SOLUTION

The preparation of the stock standard solution requires an authentic drug product for reference purposes, for example, tablets containing 50 mg of artesunate. Wrap up one reference tablet into aluminium foil and crush it down to a fine powder using a pestle. Carefully empty the aluminium foil over a 25-ml laboratory glass bottle and wash down all residual solids with 10 ml of methanol using a straight pipette. Close the bottle and shake for about three minutes until most of the solids are dissolved. Allow the solution to sit for an additional five minutes until undissolved residues settle below the supernatant liquid. The solution obtained should contain 5 mg of total artesunate per ml and be labelled as '*Artesunate Stock Standard Solution*'. Freshly prepare this solution for each test. Continue to work with the clear or hazy supernatant liquid.

IV. PREPARATION OF THE WORKING STANDARD SOLUTION 100% (UPPER WORKING LIMIT)

The artesunate stock standard solution requires no further dilution. It already represents the final working concentration of 5 mg of total drug per ml. Just for more convenient handling, some of the supernatant liquid may want to be transferred into a 10-ml vial.

This higher working standard solution represents a drug product of good quality containing 100 % of artesunate.

V. PREPARATION OF THE WORKING STANDARD SOLUTION 80% (LOWER WORKING LIMIT)

Pipette 4 ml of the stock standard solution into a 10-ml vial and add 1 ml of methanol. Close and shake the vial. The solution obtained should contain 4 mg of total drug per ml and be labelled as '*Artesunate Working Standard Solution 80%*'.

This lower working standard solution represents a drug product of poor quality containing just 80% of the amount of artesunate as stated on the product's label. In the current investigation, this drug level represents the lower acceptable limit for a given product.

VI. PREPARATION OF THE STOCK SAMPLE SOLUTION FROM A PRODUCT CLAIMING TO CONTAIN 25 MG OF ARTESUNATE PER UNIT

Take two (!) whole tablets or capsules from an appropriate drug product sampled in the field. As usual, tablets are wrapped up into aluminium foil and crushed down to a fine powder. Transfer all the powder obtained into a 25-ml laboratory glass bottle. Powder obtained from a sample capsule should be transferred directly into the bottle adding the cap and body shells last. For extraction, add 10 ml of methanol using a straight pipette, close the bottle and shake for about three minutes until most of the solids are dissolved. Allow the solution to sit for an additional five minutes until undissolved residues settle below the supernatant liquid.

50 MG OF ARTESUNATE PER UNIT

Take one whole sample tablet or capsule and extract the powder obtained with 10 ml of methanol following the procedure described above.

60 MG OF ARTESUNATE PER UNIT

Take one whole sample tablet or capsule and extract the powder obtained with 12 ml of methanol following the procedure described above.

100 MG OF ARTESUNATE PER UNIT

Take one whole sample tablet or capsule and extract the powder obtained with 20 ml of methanol following the procedure described above.

200 MG OF ARTESUNATE PER UNIT

Take one whole sample tablet or capsule and extract the powder obtained with 40 ml of methanol following the procedure described above.

All stock sample solutions produced should finally contain 5 mg of artesunate per ml and be labelled as '*Artesunate Stock Sample Solution*'. Freshly prepare these solutions for each test. Continue to work with the clear or hazy supernatant liquids.

VII. PREPARATION OF THE WORKING SAMPLE SOLUTION

Artesunate stock sample solutions require no further dilution. They already represent the final working concentration of 5 mg of artesunate per ml. If prepared from a high quality product, the sample solution should match the concentration of artesunate of the higher working standard solution produced above.

VIII. SPOTTING

Mark an origin line parallel to and about 1.5 cm from the bottom edge of the chromatoplate and apply 2 μ l of each test and standard solution prepared as shown in the picture opposite using the microcapillary pipettes supplied.

Up to five spots can be placed on a plate. Check the uniformity of all spots using UV light of 254 nm. Even if artesunate itself stays invisible, excipients and other drug compounds will show up to facilitate verification. All spots should be circular in shape and equally spaced across the origin line. Although their intensities might differ, their diameters never should. Different intensities are due to residual amounts of tablet and capsule excipients or different drug concentrations in the sample solutions. A difference in spot size, however, relates to poor spotting. Repeat this step if homogeneous spotting is not achieved first time.

IX. DEVELOPMENT

Pipette 18 ml of ethyl acetate, 4 ml of acetone and precisely 0.1 ml of glacial acetic acid into the jar being used as TLC developing chamber. Close the chamber and mix thoroughly. Line the chamber's wall with filter paper and wait for about 15 minutes thus ensuring saturation of the chamber with solvent vapour. Carefully place the loaded TLC plate into the jar. Close the jar and develop the chromatoplate until the solvent front has moved about three-quarters of the length of the plate, the developing time being about 10 minutes. Remove the plate from the chamber, mark the solvent front and allow any excess solvent to evaporate using a hot plate if necessary.

X. DETECTION

When working on fixed-dose combination medicines, it is best to check the presence of other drugs before that of artesunate. For this, expose the dried chromatoplate first to UV light of 254 nm using the battery-driven lamp supplied.

After the presence or absence of other drug compounds has been verified, the chromatoplate can be exposed to sulphuric acid staining for the detection of artesunate. For this, fill the 250-ml plastic beaker supplied with a 190 ml of methanol followed by 10 ml of concentrated sulphuric acid solution and mix gently. Allow the mix to cool down and submerge the chromatoplate into the staining solution using a pair of tweezers. Instantly remove the plate and let all surplus solution run down onto paper tissue. Wipe off residual liquid from the back of the plate and continue to dry off all staining solution on the hot plate supplied. During heating, all artesunate spots are gradually becoming visible at daylight. Use this method of detection for both, artesunate identification and quantification purposes. Note that the staining process is illustrated on page 26 of the main manual issued 2008.

After staining with sulphuric acid plus heat, a detection of other drugs in co-formulated artesunate products, for example mefloquine, is possible when subjecting the chromatoplate to UV light of 365 nm in a dark room.


XI. CHROMATOPLATE OBSERVED AT DAYLIGHT AFTER EXPOSURE TO SULPHURIC ACID AND HEAT

Run No.1:
Upper working standard
representing 100% of total
artesunate

Run No.2:
A product of good quality with
acceptable artesunate content

Run No.3:
A product of poor quality with
unacceptable low artesunate content

Run No.4:
Lower working standard
representing 80% of total
artesunate



XII. OBSERVATIONS MADE AT 254 NM BEFORE STAINING

Artesunate stays invisible and no other spots should be detected unless the sample under investigation comes as co-formulated product containing also UV visible drugs, for example sulfadoxine and sulfamethoxyprazine, both having very similar travel distances around 0.57. Other drugs arising from co-formulations (amodiaquine, mefloquine, pyronaridine, pyrimethamine) and some excipients incorporated in the different tablet or capsule formulations may cause spots emerging near or on the origin line.

XIII. OBSERVATIONS MADE AT DAYLIGHT AFTER SULPHURIC ACID STAINING

A brown spot at a travel distance of about 0.46 indicates the presence of artesunate in the test solution. No other drug from common co-formulations will show up at this stage. Any additional strong spots generated by the test solution would point to foreign drugs or artesunate degradation, the latter case being more likely when associated with a smaller principal spot. Other accessory drugs arising from co-formulations or some excipients incorporated in the different tablet or capsule formulations might cause spots emerging near or on the origin line.

XIV. RESULTS & ACTIONS TO BE TAKEN

The artesunate spot in the chromatogram obtained with the test solution must correspond in terms of colour, size, intensity, shape and travel distance to that in the chromatogram obtained with the lower and higher standard solution. This result must be obtained for each method of detection. If this is not achieved, repeat the run from scratch with a second sample. Reject the batch if the drug content cannot be verified in a third run. For a second opinion, refer additional samples to a fully-fledged drug quality control laboratory. Retain samples and put the batch on quarantine until a final decision on rejection or release has been taken. For documentation purposes, take pictures of all the readings with a digital camera turning off the flash first.

Genuine or Fake?



Fighting Counterfeit Medicines · Protecting People's Life



A charitable organisation
maintained exclusively by
Merck Darmstadt · Germany

Global Pharma Health Fund
Otto-Meißner-Straße 1
60314 Frankfurt, Germany
Phone 0049-69-962387-600
info@gphf.org · www.gphf.org



USAID
FROM THE AMERICAN PEOPLE



PROMOTING THE QUALITY OF MEDICINES

U.S. Agency for

International Development

Office of Health, Infectious Diseases
and Nutrition, Ronald Reagan Bldg.,
1300 Pennsylvania Avenue NW
Washington, DC 20523-3700, USA
Phone 001-202-712-4789
Fax 001-202-216-3702
aboni@usaid.gov · www.usaid.gov

United States Pharmacopeia

Promoting the Quality
of Medicines program
12601 Twinbrook Parkway
Rockville, MD 20852-1790, USA
Phone 001-301-816-8162
Fax 001-301-816-8374
pqm@usp.org · www.pqmusp.org